

Product Datasheet

ADORA2B Antibody (orb1247841)

Description

ADORA2B Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, IF, IHC

Immunogen

The immunogen for this antibody is: CQADVKSGNGQ

Target

ADORA2B

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: In paraffin embedded Human Kidney shows cytoplasm staining in cells of renal tubules. Recommended concentration: 2-4ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of Kelly cells. Recommended concentration: 10ug/ml.

Clonality

Polyclonal

Uniprot ID

[P29275](#)

NCBI



[NP_000667.1](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: In paraffin embedded Human Kidney shows cytoplasm staining in cells of renal tubules. Recommended concentration: 2-4ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 and U2OS cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of Kelly cells. Recommended concentration: 10ug/ml.

Expiration Date


12 months from date of receipt.



orb1247841 (2 ug/ml) staining of paraffi...



orb1247841 Immunofluorescence analysis o...



orb1247841 Immunofluorescence analysis o...