

Product Datasheet

YWHAQ Antibody (orb1247867)

Description

YWHAQ Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: C-DDRKQTIDNSQ

Target

YWHAQ

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to NP_006817.1). Recommended concentration: 0.03-0.1ug/ml.

Clonality

Polyclonal

Uniprot ID

[P27348](#)

NCBI


[NP_006817.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to NP_006817.1). Recommended concentration: 0.03-0.1ug/ml.

Expiration Date

12 months from date of receipt.



orb1247867
(0.03 ug/ml)
staining of
Huma...