



Product Datasheet

PAFAH1B1 Antibody (orb1247928)


Description

PAFAH1B1 Antibody

Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: TGSVDQTVKWECR
Target	PAFAH1B1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 40kDa band observed in Human Ovary and in Rat Ovary lysates (calculated MW of 46.6kDa according to Human NP_000421.1 and to Rat NP_113951.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cortex), Thyroid and Pancreas. Recommended concentration: 3.75ug/ml.
Clonality	Polyclonal
Uniprot ID	P43034
NCBI	NP_000421.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 40kDa band observed in Human Ovary and in Rat Ovary lysates (calculated MW of 46.6kDa according to Human NP_000421.1 and to Rat NP_113951.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cortex), Thyroid and Pancreas. Recommended concentration: 3.75ug/ml.
Expiration Date	12 months from date of receipt.



orb1247928
(3.75 ug/ml)
staining of para...



orb1247928
(3.75 ug/ml)
staining of para...



orb1247928
(0.1 ug/ml)
staining of Rat O...