





Product Datasheet

MAX Antibody (orb1247944)

Description	MAX Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC, WB	
Immunogen	The immunogen for this antibody is: C-EEPQSRKKLRMEAS	
Target	MAX	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 22kDa band observed in lysates of cell line Jurkat (calculated MW of 18.3kDa according to NP_002373.3). Recommended concentration: 0.01-0.03ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Prostate. Recommended concentration: 10ug/ml. Immunofluorescence: Expression of the protein seen in the cytoplasm and nuclei of A431 and U251 cells. Recommended concentration: 10ug/ml. Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.	orb1247944 Immunofluorescence analysis o...  orb1247944 (0.01 ug/ml) staining of Jurk...  orb1247944 (10 ug/ml) staining of paraff...
Clonality	Polyclonal	
Uniprot ID	P61244	
NCBI	NP_002373.3, NP_660087.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 22kDa band observed in lysates of cell line Jurkat (calculated MW of 18.3kDa according to NP_002373.3). Recommended concentration: 0.01-0.03ug/ml. Primary incubation 1 hour at room temperature. Immunohistochemistry: Paraffin embedded Human Prostate. Recommended concentration: 10ug/ml. Immunofluorescence: Expression of the protein seen in the cytoplasm and nuclei of A431 and U251 cells. Recommended concentration: 10ug/ml. Additional	