



Product Datasheet

GYPA Antibody (orb1247965)

Description	GYPA Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247965 (1 ug/ml) staining of Human L...
Immunogen	The immunogen for this antibody is: C-SVEIENPETSQ	
Target	GYPA	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247965 (2 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 18kDa band observed in Human Liver lysates (calculated MW of 16.4kDa according to NP_002090.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Kidney shows bright staining of erythrocytes. Recommended concentration: 2-4ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P02724	
NCBI	NP_002090.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. Western Blot: Approx 18kDa band observed in Human Liver lysates (calculated MW of 16.4kDa according to NP_002090.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: In paraffin embedded Human Kidney shows bright staining of erythrocytes. Recommended concentration: 2-4ug/ml.	
Expiration Date	12 months from date of receipt.	