



Product Datasheet

CFL2 Antibody (orb1247995)

Description	CFL2 Antibody	
Species/Host	Goat	
Reactivity	Human, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247995 (0.003 ug/ml) staining of Hum...
Immunogen	The immunogen for this antibody is: C-DIKDRSTLGEK	
Target	CFL2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	HEK293 overexpressing Human Muscle Cofil...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 18kDa band observed in Human Skeletal Muscle lysates and in Rat Heart lysates (calculated MW of 18.7kDa according to NP_068733.1 and NP_619579.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: Paraffin embedded Human Heart. Recommended concentration: 3.75ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9Y281	orb1247995 (0.001 ug/ml) staining of Rat...
NCBI	NP_068733.1 , NP_619579.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 18kDa band observed in Human Skeletal Muscle lysates and in Rat Heart lysates (calculated MW of 18.7kDa according to NP_068733.1 and NP_619579.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: Paraffin embedded Human Heart. Recommended concentration: 3.75ug/ml.	
Expiration Date	12 months from date of receipt.	