


Product Datasheet

ATF2 Antibody (orb1248017)

Description

ATF2 Antibody



Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1248017
(1 ug/ml)
staining of
HepG2 (...)

Immunogen

The immunogen for this antibody is: C-IVMAPSSSQSPSGS

Target

ATF2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75kDa observed in nuclear lysates of cell lines HeLa, HepG2 and K562, which corresponds with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at room temperature.

Clonality

Polyclonal

Uniprot ID

[P15336](#)

NCBI

[NP_001871.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75kDa observed in nuclear lysates of cell lines HeLa, HepG2 and K562, which corresponds with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at room temperature.

Expiration Date

12 months from date of receipt.