

Product Datasheet

ADAM12 Antibody (orb1248028)

Description

ADAM12 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC

Applications
Immunogen

The immunogen for this antibody is: AARPLPVSPARALC

Target

ADAM12

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes


Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: In paraffin embedded Human Prostate shows strong cytoplasm staining of a group of secretory epithelial cells. Recommended concentration: 3-5ug/ml. Paraffin embedded Human Placenta. Recommended concentration: 3.75ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm and plasma membranes of A431 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.

Clonality


Polyclonal

Uniprot ID
O43184
NCBI
NP_003465.3, NP_067673.2
Dilution Range


Peptide ELISA: antibody detection limit dilution 1:64000. Immunohistochemistry: In paraffin embedded Human Prostate shows strong cytoplasm staining of a group of secretory epithelial cells. Recommended concentration: 3-5ug/ml. Paraffin embedded Human Placenta. Recommended concentration: 3.75ug/ml. Immunofluorescence: Strong expression of the protein seen in the cytoplasm and plasma membranes of A431 cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of



orb1248028 (3.8 ug/ml) staining of paraf...



orb1248028 Immunofluorescence analysis o...



orb1248028 Flow cytometric analysis of p...