

Product Datasheet

P2RX4 antibody (PE) (orb124803)

Description	Rabbit polyclonal to P2RX4 (PE).
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human P2RX4
Target	P2RX4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124803 is the PE conjugated form of orb100036
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	Q99571
Expiration Date	12 months from date of receipt.