




Product Datasheet

dUTPase/DUT Antibody (orb1248032)

Description		dUTPase/DUT Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, IHC, WB		orb1248032 (0.1 ug/ml) staining of HEK29...
Immunogen	The immunogen for this antibody is: C-QLRFARLSEHATAPT		
Target	dUTPase / DUT		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		orb1248032 (5 ug/ml) staining of paraffi...
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 20kDa band observed in lysates of cell lines HEK293 and HeLa (calculated MW of 17.7kDa according to NP_001939.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Testis and Tonsil. Recommended concentration: 5ug/ml.		
Clonality	Polyclonal		
MW	Approx 20kDa		orb1248032 (5 ug/ml) staining of paraffi...
Uniprot ID	P33316		
NCBI	NP_001020420.1, NP_001939.1, NP_001020419.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 20kDa band observed in lysates of cell lines HEK293 and HeLa (calculated MW of 17.7kDa according to NP_001939.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Testis and Tonsil. Recommended concentration: 5ug/ml.		
Expiration Date	12 months from date of receipt.		