

Product Datasheet

BOD1L Antibody (orb1248045)

Description

BOD1L Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

The immunogen for this antibody is: C-RTRGQQRVEEAP

Target

BOD1L

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 330kDa according to NP_683692.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Clonality

Polyclonal

Uniprot ID

[Q8NFC6](#)

NCBI

[NP_683692.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 330kDa according to NP_683692.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any