


Product Datasheet

GNA12 Antibody (orb1248046)

Description

GNA12 Antibody



Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1248046
(1 ug/ml)
staining of
Human P...

Immunogen

The immunogen for this antibody is: C-KKHFPDFRGDPHR

Target

GNA12

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 37kDa band observed in Human Placenta lysates and approx 42kDa in Mouse Heart lysates (calculated MW of 37.6kDa according to Human NP_001269370.1 and 44.1kDa according to Mouse NP_034432.1). Recommended concentration: 0.3-1ug/ml.



orb1248046
(0.3 ug/ml)
staining of
Mouse...

Clonality

Polyclonal

MW

Approx 37kDa

Uniprot ID

[Q03113](#)

NCBI

[NP_001269370.1](#), [NP_001280021.1](#), [NP_001269369.1](#), [NP_031379.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 37kDa band observed in Human Placenta lysates and approx 42kDa in Mouse Heart lysates (calculated MW of 37.6kDa according to Human NP_001269370.1 and 44.1kDa according to Mouse NP_034432.1). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.