


Product Datasheet

TICAM-1/TRIF Antibody (orb1248074)

Description

TICAM-1/TRIF Antibody



orb1248074
(2 ug/ml)
staining of
Daudi l...

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: PDGATFCEDFQVP

Target

TICAM-1 / TRIF

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 70kDa band observed in lysates of cell line Daudi (calculated MW of 76.4kDa according to NP_891549.1). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

MW

Approx 70kDa

Uniprot ID

[Q8IUC6](#)

NCBI

[NP_891549.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 70kDa band observed in lysates of cell line Daudi (calculated MW of 76.4kDa according to NP_891549.1). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.