

Product Datasheet

PLA2G2A Antibody (orb1248077)

Description

PLA2G2A Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

The immunogen for this antibody is: C-SYKFSNSGSRIT

Target

PLA2G2A

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 19kDa band observed in Human Placenta lysates (calculated MW of 16.1kDa according to NP_000291.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Liver and Small Intestine. Recommended concentration: 5ug/ml.

Clonality

Polyclonal

MW

Approx 19kDa

Uniprot ID

[P14555](#)

NCBI


[NP_000291.1](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 19kDa band observed in Human Placenta lysates (calculated MW of 16.1kDa according to NP_000291.1). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Liver and Small Intestine. Recommended concentration: 5ug/ml.

Expiration Date


12 months from date of receipt.



orb1248077
(0.1 ug/ml)
staining of
Place...



orb1248077
(5 ug/ml)
staining of
paraffi...



orb1248077
(5 ug/ml)
staining of
paraffi...