

Product Datasheet

ADAM19 Antibody (orb1248262)

Description

ADAM19 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

orb1248262
(0.1 ug/ml)
staining of
Splee...

Immunogen

The immunogen for this antibody is: KIQCQSSEARPLESN

Target

ADAM19

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 125kDa band observed in Human Spleen lysates (calculated MW of 101kDa according to NP_150377.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Spleen. Recommended concentration: 5ug/ml.



orb1248262
(5 ug/ml)
staining of
paraffi...

Clonality

Polyclonal

MW

Approx 125kDa

Uniprot ID

[Q9H013](#)

NCBI

[NP_150377.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 125kDa band observed in Human Spleen lysates (calculated MW of 101kDa according to NP_150377.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Spleen. Recommended concentration: 5ug/ml.

Expiration Date

12 months from date of receipt.