

Product Datasheet

CHUK Antibody (orb1248268)

Description

CHUK Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1248268
(1 ug/ml)
staining of
HeLa nu...

Immunogen

The immunogen for this antibody is: EHDHSLSCVVPQD

Target

CHUK

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 80kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 84.6kDa according to NP_001269.3). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

MW

Approx 80kDa

Uniprot ID

[O15111](#)

NCBI

[NP_001269.3](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 80kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 84.6kDa according to NP_001269.3). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.



250kDa
150kDa
100kDa



50kDa
37kDa
25kDa
20kDa
15kDa