

Product Datasheet

apolipoprotein D Antibody (orb1248280)

Description

apolipoprotein D Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

EIA, ELISA, IHC, WB

Applications
Immunogen

The immunogen for this antibody is: EIEKIPTTFENGRC

Target

apolipoprotein D

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 25kDa band observed in Human Liver and Kidney lysates (calculated MW of 21.3kDa according to NP_001638.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Kidney and Adrenal Gland. Recommended concentration: 5ug/ml. Enzyme immunoassay: This antibody is reported to work in sandwich-type ELISA both as capturer and as reporter antibody in combination with EB12282.

Clonality





Polyclonal

MW

Approx 25kDa

Uniprot ID
[P05090](#)
NCBI
[NP_001638.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 25kDa band observed in Human Liver and Kidney lysates (calculated MW of 21.3kDa according to NP_001638.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Kidney and Adrenal Gland. Recommended concentration: 5ug/ml. Enzyme immunoassay: This antibody is reported to work in


 orb1248280
(0.3 ug/ml)
staining of
Human...

 orb1248280
(5 ug/ml)
staining of
paraffi...

 orb1248280
(5 ug/ml)
staining of
paraffi...