


Product Datasheet

ILF3 Antibody (orb1248287)

Description	ILF3 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	The immunogen for this antibody is: C-QKPSYSGSYQSHQGQ
Target	ILF3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000.Immunohistochemistry:In paraffin embedded Human Kidney shows nuclear staining in convoluted tubules. Recommended concentration, 2-4ug/ml.
Clonality	Polyclonal
Uniprot ID	Q12906
NCBI	NP_060090.2, NP_036350.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000.Immunohistochemistry:In paraffin embedded Human Kidney shows nuclear staining in convoluted tubules. Recommended concentration, 2-4ug/ml.
Expiration Date	12 months from date of receipt.



orb1248287
(2 ug/ml)
staining of
paraffi...