

Product Datasheet

SLC26A1/SAT-1 Antibody (orb1248293)

Description

SLC26A1/SAT-1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

orb1248293
(2 ug/ml)
staining of
paraffi...

Immunogen

The immunogen for this antibody is: C-QTARARHRELEATD

Target

SLC26A1 / SAT-1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:4000. Immunohistochemistry: In paraffin embedded Human Kidney shows strong basal membrane staining in convoluted tubules. Recommended concentration, 2-4ug/ml.

Clonality

Polyclonal

Uniprot ID
[Q9H2B4](#)
NCBI
[NP_071325.2](#), [NP_602297.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:4000. Immunohistochemistry: In paraffin embedded Human Kidney shows strong basal membrane staining in convoluted tubules. Recommended concentration, 2-4ug/ml.

Expiration Date

12 months from date of receipt.