



Product Datasheet

GLP1R Antibody (orb1248302)

Description	GLP1R Antibody	
Species/Host	Goat	
Reactivity	Human, Porcine, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, IF, WB	orb1248302 (1 ug/ml) staining of Human H...
Immunogen	The immunogen for this antibody is: C-QHQWDGLLSYQD	
Target	GLP1R	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248302 (1 ug/ml) staining of Rat (A)...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 55kDa band observed in Human Hippocampus and in Rat and Pig Brain lysates (calculated MW of 53.0kDa according to Human NP_002053.3, 52.9kDa according to Rat NP_036860.1 and 52.7kDa according to Pig NP_001243523.1). Recommended concentration: 0.1-1ug/ml. Primary incubation 1 hour at room temperature. Immunofluorescence: Strong expression of the protein seen in the membranes of HeLa cells. Recommended concentration: 10ug/ml. Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.	 orb1248302 Immunofluorescence analysis o...
Clonality	Polyclonal	
MW	Approx 55kDa	
Uniprot ID	P43220	
NCBI	NP_002053.3	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 55kDa band observed in Human Hippocampus and in Rat and Pig Brain lysates (calculated MW of 53.0kDa according to Human NP_002053.3, 52.9kDa according to Rat NP_036860.1 and 52.7kDa according to Pig NP_001243523.1).	