

Product Datasheet

HOXC6 Antibody (orb1248342)

Description

HOXC6 Antibody

Species/Host

Goat

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

orb1248342
(0.3 ug/ml)
staining of
Human...

Immunogen

The immunogen for this antibody is: C-PNVALNSTAYDPVRH

Target

HOXC6

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note


For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Brain (Olfactory bulb) and in Rat and Pig Brain lysates (calculated MW of 26.9kDa according to NP_004494.1). In Mouse Brain lysates the band runs like a 27kDa band. Recommended concentration: 0.03-0.1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Tonsil. Recommended concentration: 5ug/ml.



orb1248342
(0.01 ug/ml)
staining of
Mous...


Clonality

Polyclonal

MW

Approx 26kDa

Uniprot ID
P09630
NCBI
NP_004494.1
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Brain (Olfactory bulb) and in Rat and Pig Brain lysates (calculated MW of 26.9kDa according to NP_004494.1). In Mouse Brain lysates the band runs like a 27kDa band. Recommended concentration: 0.03-0.1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Tonsil. Recommended concentration: 5ug/ml.

orb1248342
(0.03 ug/ml)
staining of Pig
...

Expiration Date

12 months from date of receipt.