

Product Datasheet


CREB3L2 Antibody (orb1248352)

Description


CREB3L2 Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: DRKLSELPDGE-C
Target	CREB3L2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 27.7kDa according to NP_001240704.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows staining in select nuclei. Recommended concentration, 5-10ug/ml.
Clonality	Polyclonal
MW	Approx 26kDa
Uniprot ID	Q70SY1
NCBI	NP_001240704.1 , NP_919047.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 27.7kDa according to NP_001240704.1). Recommended concentration: 1-3ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows staining in select nuclei. Recommended concentration, 5-10ug/ml.
Expiration Date	12 months from date of receipt.

250kDa
150kDa
100kDa



orb1248352
(1 ug/ml)
staining of
Human P...



orb1248352
(5 ug/ml)
staining of
paraffi...