

Product Datasheet

CREB3 Antibody (orb1248367)

Description		CREB3 Antibody	250kDa 150kDa 100kDa A B
Species/Host	Goat		50kDa 37kDa 25kDa 20kDa 15kDa
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, IHC, WB		orb1248367 (0.3 ug/ml) staining of HepG2...
Immunogen	The immunogen for this antibody is: C-HHDHTYSLPRET		
Target	CREB3		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		