

## Product Datasheet

### Prdm9 Antibody (orb1248372)

**Description**

Prdm9 Antibody

**Species/Host**

Goat

**Reactivity**

Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

orb1248372  
(2 ug/ml)  
staining of  
Rat Ova...

**Immunogen**

The immunogen for this antibody is: C-QEHFDSQNKNDK

**Target**

Prdm9

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Form/Appearance**

Liquid

**Concentration**

500 ug/mL

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 90kDa band observed in Rat Ovary lysates (calculated MW of 92.6kDa according to NP\_001102373.2). Recommended concentration: 2-3ug/ml.

**Clonality**

Polyclonal

**MW**

Approx 90kDa

**Uniprot ID**

**P0C6Y7**

**NCBI**

**NP\_001102373.2**

**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 90kDa band observed in Rat Ovary lysates (calculated MW of 92.6kDa according to NP\_001102373.2). Recommended concentration: 2-3ug/ml.

**Expiration Date**

12 months from date of receipt.