

Product Datasheet

BEND5 Antibody (orb1248380)

Description

BEND5 Antibody

250kDa
150kDa
100kDa

A B

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1248380
(0.3 ug/ml)
staining of
Daudi...

Immunogen

The immunogen for this antibody is: CKNEERREAKYNLQ

Target

BEND5

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 48kDa band observed in Human Tonsil and Spleen and Mouse Spleen lysates and also in lysates of cell line Daudi (calculated MW of 48.2kDa according to NP_078879.2). Recommended concentration: 0.3-1ug/ml.

Clonality

Polyclonal

MW

Approx 48kDa

Uniprot ID

[Q7L4P6](#)

NCBI

[NP_078879.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 48kDa band observed in Human Tonsil and Spleen and Mouse Spleen lysates and also in lysates of cell line Daudi (calculated MW of 48.2kDa according to NP_078879.2). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.