


Product Datasheet

CIZ1 Antibody (orb1248390)

Description		CIZ1 Antibody	
Species/Host	Goat		250kDa 150kDa 100kDa 50kDa 37kDa 25kDa 20kDa 15kDa
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		orb1248390 (0.3 ug/ml) staining of HeLa ...
Immunogen	The immunogen for this antibody is: C-QPNTQDKTPSKVTAR		
Target	CIZ1		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 90kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 91.2kDa according to NP_001124490.1). Recommended concentration: 0.3-1ug/ml.		
Clonality	Polyclonal		
MW	Approx 90kDa		
Uniprot ID	Q9ULV3		
NCBI	NP_001124487.1, NP_001244905.1, NP_001244904.1, NP_036259.2, NP_001124489.1, NP_001124490.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 90kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 91.2kDa according to NP_001124490.1). Recommended concentration: 0.3-1ug/ml.		
Expiration Date	12 months from date of receipt.		