

Product Datasheet

HIC1 Antibody (orb1248399)

Description

HIC1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: NTHCAELYASGP

Target

HIC1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: In transfected HEK293 transiently expressing Human HIC1 a band of approx. 100kDa is observed. This band is not observed in the mock-transfected HEK293. The same band is observed endogenously in the cell lines A549 and Jurkat. The calculated molecular size is 74.5kDa according to NP_006488.2. Recommended concentration, 0.5-1ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q14526](#)

NCBI



[NP_001091672.1](#), [NP_006488.2](#)

Dilution Range


Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: In transfected HEK293 transiently expressing Human HIC1 a band of approx. 100kDa is observed. This band is not observed in the mock-transfected HEK293. The same band is observed endogenously in the cell lines A549 and Jurkat. The calculated molecular size is 74.5kDa according to NP_006488.2. Recommended concentration, 0.5-1ug/ml.

Expiration Date

12 months from date of receipt.



HEK293
overexpressing
Human HIC1
and pro...



orb1248399 (1
ug/ml) staining
of A549 (A...