

Product Datasheet

Pycard Antibody (orb1248434)

Description

Pycard Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

The immunogen for this antibody is: C-SVPAQSTARTGH

Target

Pycard

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Mouse Spleen, Lymph Node and Thymus lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Clonality

Polyclonal

Uniprot ID
Q9EPB4
NCBI
NP_075747.3
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Mouse Spleen, Lymph Node and Thymus lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Expiration Date

12 months from date of receipt.