

Product Datasheet

Pycard Antibody (orb1248435)

Description

Pycard Antibody

Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: C-DLTDKLVSYYLES

Target

Pycard

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 21+24kDa band observed in Mouse Spleen lysates (calculated MW of 21.5kDa according to NP_075747.3). The observed molecular weights correspond to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at room temperature.

Clonality

Polyclonal

MW

Approx 21 + 23kDa

Uniprot ID

[Q9ULZ3](#)

NCBI


[NP_075747.3](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 21+24kDa band observed in Mouse Spleen lysates (calculated MW of 21.5kDa according to NP_075747.3). The observed molecular weights correspond to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at room temperature.

Expiration Date

12 months from date of receipt.



orb1248435
(0.5 ug/ml)
staining of
Mouse...