



Product Datasheet

LAT2 Antibody (orb1248436)

Description		LAT2 Antibody		
Species/Host	Goat	 orb1248436 (0.01 ug/ml) staining of Huma...		
Reactivity	Human			
Conjugation	Unconjugated			
Tested Applications	ELISA, WB			
Immunogen	The immunogen for this antibody is: C-RSEKIYQQRSLRED			
Target	LAT2			
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.			
Form/Appearance	Liquid			
Concentration	500 ug/mL			
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.			
Note	For research use only			
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx 26kDa band observed in Human Peripheral Blood Lymphocytes lysates (calculated MW of 26.6kDa according to NP_054865.2). Recommended concentration: 0.01-0.03ug/ml.			
Clonality	Polyclonal			
MW	Approx 26kDa			
Uniprot ID	Q9GZY6			
NCBI	NP_054865.2			
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx 26kDa band observed in Human Peripheral Blood Lymphocytes lysates (calculated MW of 26.6kDa according to NP_054865.2). Recommended concentration: 0.01-0.03ug/ml.			
Expiration Date	12 months from date of receipt.			