

Product Datasheet

TRIP4/ASC1 Antibody (orb1248438)

Description

TRIP4/ASC1 Antibody

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The immunogen for this antibody is: C-DGQKSGDHLKRGR
Target	TRIP4 / ASC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Human Testis, Thyroid gland and Ovary lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Clonality	Polyclonal
Uniprot ID	Q15650
NCBI	NP_057297.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Human Testis, Thyroid gland and Ovary lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Expiration Date	12 months from date of receipt.