


Product Datasheet

ATF6 Antibody (orb1248439)

Description	ATF6 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-NSYRYDHSVSNKD	orb1248439 (2 ug/ml) staining of HeLa ly...
Target	ATF6	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 80kDa band observed in lysates of cell line HeLa (calculated MW of 74.6kDa according to NP_031374.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Rat Liver for this particular batch.	
Clonality	Polyclonal	
MW	Approx 80kDa	
Uniprot ID	P18850	
NCBI	NP_031374.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 80kDa band observed in lysates of cell line HeLa (calculated MW of 74.6kDa according to NP_031374.2). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Rat Liver for this particular batch.	
Expiration Date	12 months from date of receipt.	