

Product Datasheet

BACH1 Antibody (orb1248440)

Description	BACH1 Antibody	 250kDa 150kDa 100kDa 50kDa 37kDa 25kDa 20kDa 15kDa
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1248440 (1 ug/ml) staining of Human P...
Immunogen	The immunogen for this antibody is: HASVQPNERSENEC	
Target	BACH1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 80kDa band observed in Human Spleen and Peripheral Blood Lymphocytes lysates (calculated MW of 82.0kDa according to NP_001177.1). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
MW	Approx 80kDa	
Uniprot ID	O14867	
NCBI	NP_001177.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 80kDa band observed in Human Spleen and Peripheral Blood Lymphocytes lysates (calculated MW of 82.0kDa according to NP_001177.1). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	