


Product Datasheet

CD19 Antibody (orb1248441)

Description		CD19 Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		orb1248441 (2 ug/ml) staining of Jurkat ...
Immunogen	The immunogen for this antibody is: C-ENDSNLGQDQLS		
Target	CD19		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 61.1kDa according to NP_001761.3). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3ug/ml.		
Clonality	Polyclonal		
MW	Approx 75kDa		
Uniprot ID	P15391		
NCBI	NP_001761.3 , NP_001171569.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 61.1kDa according to NP_001761.3). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3ug/ml.		
Expiration Date	12 months from date of receipt.		