

Product Datasheet

CSTF3 Antibody (orb1248477)

Description	CSTF3 Antibody	 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa
Reactivity	Human	 25kDa 20kDa
Conjugation	Unconjugated	 15kDa
Tested Applications	ELISA, WB	orb1248477 (0.1 ug/ml) staining of nucle...
Immunogen	The immunogen for this antibody is: C-RETKGKLPSYKE	
Target	CSTF3	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 75kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 82.9kDa according to NP_001317.1). The observed molecular weight corresponds to the expected size of 77kDa. Recommended concentration: 0.03-0.1ug/ml.	
Clonality	Polyclonal	
MW	Approx 75kDa	
Uniprot ID	Q12996	
NCBI	NP_001317.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 75kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 82.9kDa according to NP_001317.1). The observed molecular weight corresponds to the expected size of 77kDa. Recommended concentration: 0.03-0.1ug/ml.	
Expiration Date	12 months from date of receipt.	