

Product Datasheet

**Serotonin Transporter antibody (PE) [Out of stock]
(orb124850)**

Description

Rabbit polyclonal antibody against

PE

Species/Host

Rabbit

Conjugation

PE

Immunogen

KLH conjugated synthetic peptide derived from human 5-HTT

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Application notes

orb124850 is the PE conjugated form of orb13226

Isotype

IgG

Clonality

Polyclonal

MW

70kDa

Uniprot ID
P31645, Q60857
Expiration Date

12 months from date of receipt.