

## Product Datasheet

### ADRA1A Antibody (orb1248509)

**Description**

ADRA1A Antibody

**Species/Host**

Goat

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

orb1248509  
(1 ug/ml)  
staining of  
Human P...

**Immunogen**

The immunogen for this antibody is: C-PSLDKNHQVPTIK

**Target**

ADRA1A

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Form/Appearance**

Liquid

**Concentration**

500 ug/mL

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in Human Heart and Prostate lysates (calculated MW of 51.5kDa according to NP\_000671.2). Recommended concentration: 1-3ug/ml.

**Clonality**

Polyclonal

**MW**

Approx 55kDa

**Uniprot ID**
**P35348**
**NCBI**
**NP\_000671.2**
**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in Human Heart and Prostate lysates (calculated MW of 51.5kDa according to NP\_000671.2). Recommended concentration: 1-3ug/ml.

**Expiration Date**

12 months from date of receipt.