



Product Datasheet

OSMR Antibody (orb1248546)

Description		OSMR Antibody	
Species/Host	Goat		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, FC, WB		orb1248546 (1 ug/ml) staining of U937 (A...
Immunogen	The immunogen for this antibody is: C-QIHGEQLDPHVT		
Target	OSMR		
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.		
Form/Appearance	Liquid		
Concentration	500 ug/mL		
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.		orb1248546 Flow cytometric analysis of p...
Note	For research use only		
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 125kDa band observed in lysates of cell line U937 and approx 70+125kDa in lysates of cell line A431 and HepG2 (calculated MW of 111kDa according to NP_003990.1 (isoform 1) and NP_001310435.1 (isoform 3), and 66.0kDa according to NP_001310436.1 (isoform 4)). Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at room temperature. Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.		
Clonality	Polyclonal		
MW	Approx 125kDa		
Uniprot ID	Q99650		
NCBI	NP_003990.1, NP_001161827.1, NP_001310436.1, NP_001310435.1		
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 125kDa band observed in lysates of cell line U937 and approx 70+125kDa in lysates of cell line A431 and HepG2 (calculated MW of 111kDa according to NP_003990.1 (isoform 1) and NP_001310435.1 (isoform 3), and 66.0kDa according to NP_001310436.1 (isoform 4)). Recommended concentration: 0.3-1ug/ml. Primary incubation 1 hour at		

orb1248546
(1 ug/ml)
staining of
U937 (A...



orb1248546
Flow
cytometric
analysis of
p...