


Product Datasheet

TSPAN14 Antibody (orb1248563)

Description	TSPAN14 Antibody	
Species/Host	Goat	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1248563 (2 ug/ml) staining of Mouse L...
Immunogen	The immunogen for this antibody is: C-SDIEAVKAGHH	
Target	TSPAN14	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 30kDa band observed in Mouse Lung lysates (calculated MW of 30.6kDa according to Mouse NP_666040.1). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
MW	Approx 30kDa	
Uniprot ID	Q8NG11	
NCBI	NP_112189.2 , NP_001121781.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 30kDa band observed in Mouse Lung lysates (calculated MW of 30.6kDa according to Mouse NP_666040.1). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	