


Product Datasheet

SLC26A6 Antibody (orb1248570)

Description		 250kDa 150kDa 100kDa
Species/Host	Goat	 50kDa 37kDa
Reactivity	Human	 25kDa 20kDa
Conjugation	Unconjugated	 15kDa
Tested Applications	ELISA, IHC, WB	orb1248570 (0.1 mg/ml) overnight stainin...
Immunogen	The immunogen for this antibody is: C-SAPRTHQWRT	
Target	SLC26A6	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	 250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248570 (10 ug/ml) staining of Human ...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 90kDa band observed in Human renal cortical brush-border membrane lysates (calculated MW of 83.0kDa according to NP_075062.2). Recommended concentration: 10ug/ml for overnight incubation. Immunohistochemistry: In paraffin embedded Human Kidney shows staining of brush borders along the entire proximal convoluted tubule. Recommended concentration, 0.1-0.3mg/ml. Please note the unusually low dilution required for this type of application!	
Clonality	Polyclonal	
MW	Approx 90kDa	
Uniprot ID	Q9BX59	
NCBI	NP_599025.2 , NP_075062.2 , NP_001035544.1 , NP_602298.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 90kDa band observed in Human renal cortical brush-border membrane lysates (calculated MW of 83.0kDa according to NP_075062.2). Recommended concentration: 10ug/ml for overnight incubation. Immunohistochemistry: In paraffin embedded Human Kidney shows staining of brush borders along the entire proximal convoluted tubule. Recommended concentration, 0.1-0.3mg/ml. Please note the unusually low dilution required for this type of application!	