


Product Datasheet

BDNF Antibody (orb1248588)

Description

BDNF Antibody


Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	The immunogen for this antibody is: ETKCNPMGYTKE
Target	BDNF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:2000.Immunofluorescence: Expression of the protein seen in the nucleus and cytoplasm of Glioblastoma U251 cells. Recommended concentration: 10ug/ml.
Clonality	Polyclonal
Uniprot ID	P23560
NCBI	NP_733931.1 , NP_733928.1 , NP_733929.1 , NP_001700.2 , NP_733930.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000.Immunofluorescence: Expression of the protein seen in the nucleus and cytoplasm of Glioblastoma U251 cells. Recommended concentration: 10ug/ml.
Expiration Date	12 months from date of receipt.



orb1248588 (1 ug/ml) staining of PFA-per...



orb1248588 Immunofluorescence analysis o...



orb1248588 (5 ug/ml) staining of paraffi...