



Product Datasheet

ANTXR2/CMG2 Antibody (orb1248590)

Description	ANTXR2/CMG2 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	The immunogen for this antibody is: C-HEGLKLANEQIQK
Target	ANTXR2 / CMG2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000.Immunohistochemistry:Paraffin embedded Human Liver and Kidney. Recommended concentration: 2.5ug/ml.
Clonality	Polyclonal
MW	Approx 60kDa
Uniprot ID	P58335
NCBI	NP_001139266.1, NP_477520.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000.Immunohistochemistry:Paraffin embedded Human Liver and Kidney. Recommended concentration: 2.5ug/ml.
Expiration Date	12 months from date of receipt.



orb1248590
(2.5 ug/ml)
staining of
paraf...



orb1248590
(2.5 ug/ml)
staining of
paraf...