



Product Datasheet

CHRNA7 Antibody (orb1248594)

Description	CHRNA7 Antibody	
Species/Host	Goat	
Reactivity	Human, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	
Immunogen	The immunogen for this antibody is: KRPGEDKVRPACQHKQ	orb1248594
Target	CHRNA7	(5 ug/ml)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	staining of paraffi...
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Immunohistochemistry: In paraffin embedded Rat Cerebral Cortex and Rat Hippocampus shows membranous staining in neural cells. Recommended concentration: 5-10ug/ml. Paraffin embedded Human Brain (Cortex). Recommended concentration: 3-6ug/ml.	 orb1248594 (3 ug/ml) staining of paraffi...
Clonality	Polyclonal	
Uniprot ID	P36544	
NCBI	NP_000737.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Immunohistochemistry: In paraffin embedded Rat Cerebral Cortex and Rat Hippocampus shows membranous staining in neural cells. Recommended concentration: 5-10ug/ml. Paraffin embedded Human Brain (Cortex). Recommended concentration: 3-6ug/ml.	
Expiration Date	12 months from date of receipt.	