


Product Datasheet

GPI/Neuroleukin Antibody (orb1248598)


Description

GPI/Neuroleukin Antibody


Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: C-RERMFNGEKINYT
Target	GPI / Neuroleukin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 60kDa band observed in Human Heart lysates and approx 55kDa observed in lysates of cell lines HeLa, A431, A549, MCF7, Jurkat and K562 (calculated MW of 63.1kDa according to NP_000166.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Testis, Thyroid and Kidney. Recommended concentration: 3.75ug/ml.
Clonality	Polyclonal
MW	Approx 60kDa
Uniprot ID	P06744
NCBI	NP_000166.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 60kDa band observed in Human Heart lysates and approx 55kDa observed in lysates of cell lines HeLa, A431, A549, MCF7, Jurkat and K562 (calculated MW of 63.1kDa according to NP_000166.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Testis, Thyroid and Kidney. Recommended concentration: 3.75ug/ml.



orb1248598
(0.3 ug/ml)
staining of
Human...



orb1248598
(0.05 ug/ml)
staining of
HeLa...



orb1248598
(3.75 ug/ml)
staining of
para...