

Product Datasheet

STIM1 Antibody (orb1248617)

Description

STIM1 Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: C-YIQNRYSKEHMKK
Target	STIM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Successful use of this antibody was reported in Western blot showing an approx. 80kDa in Human platelets: PMID 18439569
Clonality	Polyclonal
MW	Approx 80kDa
Uniprot ID	Q13586
NCBI	NP_003147.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Successful use of this antibody was reported in Western blot showing an approx. 80kDa in Human platelets: PMID 18439569
Expiration Date	12 months from date of receipt.