


Product Datasheet

EHD2 Antibody (orb1248621)

Description	EHD2 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, IF, WB	orb1248621 (0.1 µg/ml) staining of A549...
Immunogen	The immunogen for this antibody is: CRLVPPSKRRHKGSA	
Target	EHD2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248621 (0.3 µg/ml) staining of Huma...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Human Placenta lysates (calculated MW of 61.2kDa according to NP_055416.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunofluorescence: strong expression of the protein seen in the plasma membrane of A431 cells and the nuclei of HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9NZN4	
NCBI	NP_055416.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Human Placenta lysates (calculated MW of 61.2kDa according to NP_055416.2). Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunofluorescence: strong expression of the protein seen in the plasma membrane of A431 cells and the nuclei of HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.	
Expiration Date	12 months from date of receipt.	