

Product Datasheet

KLF2 / LKLF Antibody [Out of stock] (orb1248622)

Description

KLF2 / LKLF Antibody [Out of stock]

250kDa
150kDa
100kDa

Species/Host

Goat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The immunogen for this antibody is: ALSEPILPSFST-C

Target

KLF2 / LKLF

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 37kDa band observed in Human Heart and Lung lysates (calculated MW of 37.4kDa according to NP_057354.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Rat Brain and Mouse Heart for this particular batch.

Clonality

Polyclonal

MW

Approx 37kDa

Uniprot ID

[Q9Y5W3](#)

NCBI

[NP_057354.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 37kDa band observed in Human Heart and Lung lysates (calculated MW of 37.4kDa according to NP_057354.1). Recommended concentration: 1-3ug/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Rat Brain and Mouse Heart for this particular batch.

Expiration Date

12 months from date of receipt.



orb1248622
(1 ug/ml)
staining of
Human H...