


Product Datasheet

SNX8 Antibody (orb1248623)

Description	SNX8 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	
Immunogen	The immunogen for this antibody is: TPPCSPPEDGLCPH	orb1248623 (4 ug/ml) staining of paraffi...
Target	SNX8	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000.Immunohistochemistry:In paraffin embedded Human Kidney shows vesicular staining in the cytoplasm of preferentially DCT epithelial cells. Recommended concentration, 4-8ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9Y5X2	
NCBI	NP_037453.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000.Immunohistochemistry:In paraffin embedded Human Kidney shows vesicular staining in the cytoplasm of preferentially DCT epithelial cells. Recommended concentration, 4-8ug/ml.	
Expiration Date	12 months from date of receipt.	