


Product Datasheet

NUMBL Antibody (orb1248625)

Description

NUMBL Antibody

Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: SRSAAASGGPRRPE-C
Target	NUMBL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: In transfected HEK293 transiently expressing Human Numbl and DYKDDDDK tagged), a band of approx. 85kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration, 0.5-1ug/ml.
Clonality	Polyclonal
Uniprot ID	Q9Y6R0
NCBI	NP_004747.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: In transfected HEK293 transiently expressing Human Numbl and DYKDDDDK tagged), a band of approx. 85kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration, 0.5-1ug/ml.
Expiration Date	12 months from date of receipt.



HEK293
lysate (10 ug
protein in
RIPA buf...