


Product Datasheet

HTR2A Antibody (orb1248628)

Description	HTR2A Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	
Immunogen	The immunogen for this antibody is: C-QRSIHREPGSYTGR	orb1248628 (2 ug/ml) staining of paraffi...
Target	HTR2A	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248628
Note	For research use only	Negative Control showing stai...
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Human Brain lysates of several regions gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 2ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P28223	
NCBI	NP_001159419.1 , NP_000612.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Preliminary experiments in Human Brain lysates of several regions gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 2ug/ml.	
Expiration Date	12 months from date of receipt.	